

# Panther TLSL 464



## Quad Control TLSO

The Panther TLSL spinal brace delivers maximum spinal control. The posterior-lateral shell, rigid anterior plate, and twin lateral plates overlap to provide circumferential support & control.

With its choice of pectoral pads or sternal plate and attached shoulder straps, the Panther TLSL offers triplanar control from the symphysis pubis to the sternum anteriorly and from the gluteal fold to the scapular spines posteriorly.

The string-and-pulley compression system allows for easy cinching, even by elderly arthritic hands.

The Panther TLSL is universal. It requires no tools to fit. There's no need to cut or fold a belt. Just add or remove belt segments.

PDAC Approved: L0464



SCAN FOR  
PANTHER CATALOG



The Best Value in Spinal Bracing

800-867-6070 | [www.bluedortho.com](http://www.bluedortho.com)

# Panther TLSL 464



## SPECIFICATIONS

- Universal Panther belt
- Lateral overlapping plastic
- Smooth dual-pull compression
- Sacral to sternal spine control
- Ant. & Post. height adjustments
- Shoulder/axilla straps

## UNIVERSAL SIZE

- One size fits all.
- 25" to 66" and beyond! Just add extra belt segments.

## PDAC

- L0464 Custom Fit only

## LORDOSIS

The posterior polymer plate is preset at 15°. It may be heat molded for exact fit during custom fitting.

STERNAL PLATE TOP



TLSL-464-I

PECTORAL PAD TOP



TLSL-464-FP

