

LEOPARD ONE TLSL



One means one size!

The Leopard TLSL spinal brace delivers maximum spinal control. The posterior, anterior and long bi-lateral control plates overlap to provide circumferential support & control.

With its sternal STEP and attached shoulder straps, the Leopard TLSL offers triplanar control from the symphysis pubis to the sternal notch anteriorly and from the gluteal fold to the scapular spines posteriorly.

The quad-strap compression system snugs the brace firmly.

To fit the Leopard TLSL, set the torso circumference by trimming and reattaching the lateral extensions, then replace the straps. Adjust the long lateral pouch and the shoulder straps, then don the brace, snug it to the torso, and set the STEP height.

PDAC Approved: L0464



The Best Value in Spinal Bracing

800-867-6070 | www.bluedortho.com

LEOPARD ONE TLSL



SPECIFICATIONS

- Quad Control TLSO
- Posterior-lateral shell
- Plush tri-laminate fabric
- Smooth quad-pull compression
- Sacral to thoracic spinal control
- Adjustable shoulder/axilla straps
- Choice of STEP
 - I – Greater vertical range
 - II – Quick-adjust extension

PDAC APPROVED

- L0464 Custom Fit

PRECISION ADJUSTMENT

The Leopard TLSL anterior plate, posterior-lateral shell, lateral control plates, and sternal STEP plate can be reformed with a heat gun for a precise fit.

UNIVERSAL SIZE

One size fits 30" to 65"

