

LEOPARD ONE TL



One means one size!

The Leopard TL spinal brace is an excellent TLSO for the treatment of compression fractures, excessive kyphosis, or for post-fusion stabilization.

The quad-strap compression system snugs the brace firmly.

The Leopard TL is easy to fit. Set the size by trimming and reattaching the lateral extensions, reset the straps, adjust the shoulder straps, don and tighten.

That's it!

PDAC Approved: L0456/L0457



The Best Value in Spinal Bracing

800-867-6070 | www.bluedortho.com

LEOPARD ONE TL



SPECIFICATIONS

- Extension Assist TLSO
- Plush tri-laminate fabric
- One size – easy to stock
- Conforms to any body shape
- Adjustable-height posterior strut
- Adjustable shoulder/axilla straps
- Sacral to thoracic flexion control
- Adjustable quad-pull compression

UNIVERSAL SIZE

One size fits 30" to 65"

PDAC APPROVED

- L0456 Custom Fit
- L0457 Off the Shelf (OTS)

PRECISION ADJUSTMENT

The Leopard TL's genuine Kydex anterior and posterior control plates can be reformed with a heat gun for a precise fit.

