

# Lynx LSO



The Lynx spinal brace is made of reinforced elastic with polymer plates in pockets. Two-stage elastic cinching allows for easy compression. You might say, "It's a cinch to cinch."

The orthopedic latex-free elastic offers ventilation for cool comfort as it contours itself to any body shape.

The molded polymer control plates are heat remoldable.

The Lynx MP (pictured this page) is available in two heights, standard and tall.

The Lynx SP (pictured on back) is identical to the standard Lynx MP, except it adds a tall LSO back insert.

Lynx MP - PDAC approved L0627/L0642

Lynx SP - PDAC approved L0631/L0648



Lynx MP



Lynx MP



The Best Value in Spinal Bracing

800-867-6070 | [www.bluedortho.com](http://www.bluedortho.com)

# Lynx LSO



## SPECIFICATIONS

- Double-thickness breathable elastic
- Conforms to any body shape
- Dual cinch straps with pull handles

## SIZES

- Small 26" - 32"
- Medium 30" - 37"
- Large 34" - 41"
- X-Large 38" - 47"
- 2X-Large 44" - 54"
- 3X-Large 50" - 62"
- 4X-Large 56" - 68"

## PDAC APPROVED

- Lynx MP: L0627 Custom Fit
- Lynx MP: L0642 Off the Shelf (OTS)
- Lynx SP: L0631 Custom Fit
- Lynx SP: L0648 Off the Shelf (OTS)

## OPTIONS

- 0°, 5°, **15°**, or 25° lordosis
- **Standard** or Pendulous anterior  
**Bolded** is default configuration



**LYNX SP**

