

Puma Una SL



Two Sizes fit all!

The Puma Una SL delivers comfort and control for the patient with moderate to severe low-back pain. The plush breathable fabric provides the comfort. The rigid Kydex provides the control.

The dual cinch straps allow the brace to be smoothly snugged, even by elderly patients with diminished hand strength.

With just two sizes, Petite and Grande, the Puma Una spinal brace can accurately fit torsos from 30" to 66".

PDAC Approved: L0637/L0650



The Best Value in Spinal Bracing

800-867-6070 | www.bluedortho.com

Puma Una SL



SPECIFICATIONS

- Easy-pull strap closure
- Provides support from S1 to T9
- Lateral plates for coronal control
- Tall back plate for sagittal control
- Plush breathable tri-laminate fabric

SIZES

Two sizes fit all

- Petite - 30" to 44"
- Grande - 42" to 66"

PDAC APPROVED

- L0637 Custom Fit
- L0650 Off the Shelf (OTS)

OPTIONS

- 0°, 5°, **15°**, or 25° lordosis
- **Standard** or Pendulous anterior
- **Standard** or STEP-ready (thicker, pre-drilled) front plate

***Bolded** is default configuration*

