



Puma Una Spinal Brace





Dual Coding

Most BDO spinal braces are dual-coded by PDAC. That means that are designed to be dispensed as is (OTS) or to be custom fitted by an experienced professional.

The polymer plates may be heat remolded to conform to the torso. For other custom modifications, see www.bluedortho.com.

Satisfaction

When you fit a Blue Diamond Orthopedic brace, you'll be rewarded with the gratitude of trusting patients and the respect of referring physicians.

Ask Us About

- Custom-branding
- Becoming a sales representative

Warranty

Blue Diamond Orthopedic warrants its spinal braces for one year.

Puma Una MP



Two Sizes fit all!

The Puma Una MP delivers comfort and control for the patient with mild to moderate low-back pain. The plush breathable fabric provides the comfort. The rigid Kydex provides the control.

The dual cinch straps allow the brace to be smoothly snugged, even by elderly patients with diminished hand strength.

With just two sizes, Petite and Grande, the Puma Una spinal brace can accurately fit torsos from 30" to 66".

PDAC Approved: L0627/L0642



The Best Value in Spinal Bracing

800-867-6070 | www.bluedortho.com

Puma Una MP



SPECIFICATIONS

- Easy-pull strap closure
- Short back plate for mild control
- Provides support from sacrum to L1
- Plush breathable tri-laminate fabric

SIZES

Two sizes fit all

- Petite - 30" to 44"
- Grande - 42" to 66"

PDAC APPROVED

- L627 Custom Fit
- L642 Off the Shelf (OTS)

OPTIONS

- 0°, 5°, **15°**, or 25° lordosis
- **Standard** or Pendulous anterior
- **Standard** or STEP-ready (thicker, pre-drilled) front plate

***Bolded** is default configuration*



Puma Una SP



Two Sizes fit all!

The Puma Una SP delivers comfort and control for the patient with moderate low-back pain. The plush breathable fabric provides the comfort. The rigid Kydex provides the control.

The dual cinch straps allow the brace to be smoothly snugged, even by elderly patients with diminished hand strength.

With just two sizes, Petite and Grande, the Puma Una spinal brace can accurately fit torsos from 30" to 66".

PDAC Approved: L0631/L0648



The Best Value in Spinal Bracing

800-867-6070 | www.bluedortho.com

Puma Una SP



SPECIFICATIONS

- Easy-pull strap closure
- Provides support from S1 to T9
- Tall back plate for sagittal control
- Plush breathable tri-laminate fabric

SIZES

Two sizes fit all

- Petite - 30" to 44"
- Grande - 42" to 66"

PDAC APPROVED

- L0631 Custom Fit
- L0648 Off the Shelf (OTS)

OPTIONS

- 0°, 5°, **15°**, or 25° lordosis
- **Standard** or Pendulous anterior
- **Standard** or STEP-ready (thicker, pre-drilled) front plate

***Bolded** is default configuration*



Puma Una SL



Two Sizes fit all!

The Puma Una SL delivers comfort and control for the patient with moderate to severe low-back pain. The plush breathable fabric provides the comfort. The rigid Kydex provides the control.

The dual cinch straps allow the brace to be smoothly snugged, even by elderly patients with diminished hand strength.

With just two sizes, Petite and Grande, the Puma Una spinal brace can accurately fit torsos from 30" to 66".

PDAC Approved: L0637/L0650



The Best Value in Spinal Bracing

800-867-6070 | www.bluedortho.com

Puma Una SL



SPECIFICATIONS

- Easy-pull strap closure
- Provides support from S1 to T9
- Lateral plates for coronal control
- Tall back plate for sagittal control
- Plush breathable tri-laminate fabric

SIZES

Two sizes fit all

- Petite - 30" to 44"
- Grande - 42" to 66"

PDAC APPROVED

- L0637 Custom Fit
- L0650 Off the Shelf (OTS)

OPTIONS

- 0°, 5°, **15°**, or 25° lordosis
- **Standard** or Pendulous anterior
- **Standard** or STEP-ready (thicker, pre-drilled) front plate

***Bolded** is default configuration*

